

## 4. Digestions

### 4.1 Microwave Digestions

4.1.1 Estimation of Weight

Dana Zimmer, Rhena Schumann

### Question:

It has to be estimated, which weight of sample material is necessary and in which range the P concentration is in the extract. Th estimation results from estimated P concentrations in environmental samples (see chapter 1) and a standard digestion method.

### Known requirements:

- An estimated P concentration range of the sample has to be known (see chapter 1): e.g. in bone char 100...150 g P kg<sup>-1</sup>. In the following formula it is calculated with 100 g P to make it easier. Calculation for 150 g P are analogue.
- A standard procedure including supposed weigh-in, dilution and so on is selected. Bone char has an organic matrix and can therefore be digested according to a standard method for plant material. Normally, plant material is processed as in the following:
  - Weigh in 0.1 g sample
  - Digest with 5 ml conc. HNO<sub>3</sub> and 3 ml 30 %  $H_2O_2$
  - Fill to 100 ml with ultra-pure water
  - Measure P at ICP-OES
  - ► The middle standard has a concentration of 10 mg P per litre
  - the P concentration in the extract should be in the range of the middle standard of calibration line and/or not exceed highest standard by a factor of 10



# Stepwise procedure to estimate P concentration in the extract according to the standard procedure:

the estimated P concentration in the sample (100...150 g P per 1000 g) is converted to the standard weigh-in by a ratio equation:

$$\frac{100 \ g \ P}{x} \equiv \frac{1000 \ g}{0.1 \ g}$$

conversion results in the following formula:

$$\frac{100 \ g \ P \times 0.1 \ g \ weight}{1000 \ g} = x = \ 0.01 \ g \ P$$

In 0.1 g bone char there are 0.01 to 0.015 g P.

- This amount of P in the weigh-in mass corresponds to the P amount in the extract after microwave digestion (from this weigh-in). If after the microwave digestion the extract with HNO<sub>3</sub> and H<sub>2</sub>O<sub>2</sub> is filled to 100 ml with ultra-pure water, between 0.01 to 0.015 g P can be expected in the 100 ml extract.
- This P amount in 100 ml extract is converted by the ratio equation to P concentrations per litre:

$$\frac{0.01 \ g \ P}{x} \equiv \frac{100 \ ml}{1000 \ ml}$$

conversion results in the following formula:

 $\frac{0.01 \ g \ P \times 1000 \ ml}{100 \ ml} = x = \ 0.1 \ g \ P \ per \ litre$ 

► The extract has a P concentration between 0.1 and 0.15 g P per litre. This corresponds to a concentration of 100 to 150 mg P per litre in the extract.

### **Comparison of the estimated P concentration to the standards:**

- ► Comparison
  - ▶ P concentration of the middle standard: 10 mg P per litre
  - ▶ P concentration of the extract 100 to 150 mg P per litre



- For a standard weigh-in of 0,1 g bone char and filling to 100 ml the P concentration in the extract is 10 to 15 times as high as the middle standard.
- Conclusion
  - ▶ P concentration in the extract should be by a fifth to tenth lower!
- ► Generally, there are 2 opportunities to achieve this:
  - Decreasing the weigh-in to a tenth: weigh-in only 0.01 g bone char instead of 0.1 g
  - Dilution of extracts by a factor of 5 to 10

### **Evaluation of pros and cons of both opportunities**

Lower weigh-in:

- Advantages
  - Less sample material is necessary.
  - No further dilution of extracts is necessary, which might cause dilution errors.
  - Less extraction agents are necessary.
- Disadvantages
  - If the material can statically charge, lower weigh-in can be problematic.
  - For heterogenic material lower weigh-in can increase standard deviation. Either the material has to be homogenised (e.g. milling) or the number of weight-ins has to be increased.
  - If other elements than P shall be measured, their concentration could fall below the detection limit if the weigh-in decreased.

#### 5- to 10-times dilution

- Advantages
  - No/less problems, which can be caused by lower weigh-in (statically charge, heterogeneity)
  - Other necessary elements are in the measurement range or by different dilution the ideal concentration ranges of different elements can be achieved.
- Disadvantages
  - Dilution error
  - ► Higher amounts of chemicals are needed.



**Decision with weighting of individual points** (most important first)

- A lower weigh-in should be considered, if: 1. the material is relatively homogenic, 2. the material is not statically charged during weighingin, 3. No other elements have to be determined in the extract and 4. less sample material is available.
- A dilution by factor 5 to 10 should be selected, if: 1. other elements have to be determined, 2. the material is heterogenic, 3. static charges could cause problems during weighing-in and 4. sufficient material is available.

### Suggestions for adjustment of the digestion method for bone char

The decrease of weigh-in by a factor of ten decreases, at a standard volume of 100 ml, the P concentration in the extract (from 100 to 150 mg  $l^{-1}$ ) to 10 to 15 mg P per litre and less chemicals are necessary for the extract (Tab. 4.1.1-1). If the volume is filled not to 100 but to 50 ml, the P concentration increased to 20 to 30 mg P per litre. This is still in the measurement range for the ICP-OES, and enables simultaneously to determine other elements. If trace elements such as Cd, Cu and Zn have to be determined, the end-volume can be decreased to 20 to 25 ml. Under theses circumstances the P concentration increased to 40 to 60 mg P per litre. That means that a further dilution of the extracts could be necessary for P measurement.

method	originally	adjusted
weigh-in	0.1 g	0.01 g
extraction	5 ml conc. HNO $_3$ and 3 ml	2.5 ml conc. HNO <sub>3</sub> and
chemicals	H <sub>2</sub> O <sub>2</sub>	1.5 ml H <sub>2</sub> O <sub>2</sub>
end volume	100 ml	50 ml
expected P concentration	100 to 150 mg P per litre	20 to 30 mg P per litre

Table 4.1.1-1 Comparison of methods for extraction of bone	char with the original and
the adjusted method	

Handbook on the selection of methods for digestion and determination of total P in environmental samples

**For citation:** Zimmer D, Baumann K, Berthold M, Schumann R (*insert year of download*): Handbook on the Selection of Methods for Digestion and Determination of Total Phosphorus in Environmental Samples. Version 1.0. DOI: 10.12754/misc-2020-0001